

**NEW!**  
**ASTM Standard Test  
Method D8148  
for Haze in Fuels!**

# Color Choice<sub>h3</sub>

- Automatic Color, Haze and Clarity measurement
- ASTM D8148 Haze in Fuels
- Eliminates operator bias
- Performs the most widely used Color and Haze scales in the Petroleum industry in just seconds
- Color tests performed ASTM D1500, Saybolt D156, D6045, D1209, D5386, ISO 2049, IP 196
- Correlates to ASTM D 4176 Procedure 2 Haze rating
- Instrument Haze Rating (IHR) 1-6
- Proprietary Haze Clarity Index (HCI) 50-100 Quantification
- Haze Trending-Time/Temp/HCI/IHR
- Test by Temperature
- State of the Art Spectrophotometer and Optics
- No moving parts or filters to degrade
- On site calibration
  - No internet connection required
- UL, TUV, CE and CSA Certified
- Simple touch screen.
  - Option for keyboard and mouse
- Store and print hundreds of results
- Easy to use, calibrate and maintain
- Output to LIMS, printer, PC or external storage device
- Proudly made in the USA



Available exclusively from:



**Choice Analytical, Inc.**

*Support You Can Trust*

For product and sales information contact: [sales@choiceanalytical.com](mailto:sales@choiceanalytical.com) or 832-271-7448  
[www.choiceanalytical.com](http://www.choiceanalytical.com)