

**NEW!**  
**ASTM Standard Test  
Method D8148  
for Haze in Fuels!**

*Clarity*  
**Choice**  
*h3*

- Automatic Haze and Clarity measurements
- ASTM D8148 Haze in Fuels
- Eliminates operator bias
- Correlates to ASTM D4176 Procedure 2 Haze rating
- Instrument Haze Rating (IHR) 1-6
- Proprietary Haze Clarity Index (HCI) 50-100 Quantification
- Haze Trending-Time/Temp/HCI/IHR
- Test by Temperature
- On site calibration
  - *No internet connection required*
- State of the Art Optics
- No moving parts or filters to degrade
- UL, TUV, CE and CSA Certified
- Simple touch screen.
  - Option for keyboard and mouse*
- Store and print hundreds of results
- Easy to use, calibrate and maintain
- Output to LIMS, printer, PC or external storage device
- Proudly made in the USA



*Available exclusively from:*



**Choice Analytical, Inc.**

*Support You Can Trust*

For product and sales information contact: [sales@choiceanalytical.com](mailto:sales@choiceanalytical.com) or 832-271-7448  
[www.choiceanalytical.com](http://www.choiceanalytical.com)