

**Maual Test Apparatus for Copper Strip
 Corossion Test - ASTM D130 | ASTM D4048**
 Stainless Steel Liquid (Oil) Bath for both Pressure
 Vessel & test tube Procedures.

Model: CSC-D130-L-8



Corossion Series



**Maual Test Apparatus for Silver Strip Corrosion
 Test - IP227 and Related Methods.**
 Stainless Steel Liquid (Oil) Bath
Model: SLVC-IP227-L-6


Conform to ASTM D130 | ASTM D4048

Features:

- Stainless steel bath
- Temperature : Min 35'C & Max 110'c
- Bath Volume : Approx : From 18 litre
- Bath Stability of $\pm 0.5^{\circ}\text{C}$
- Digital Temperature controller with 0.1°C resolution
- A float switch for low bath fluid level safety cut-off.
- Lid cover with holes, for test tubes and pressure vessels

Accessories Available :

- Corrosion pressure vessels
- Silicon carbide paper - 240/150 grit
- Copper and Silver test strips
- Glass test tubes for D130 & IP227
- ASTM corrosion code standard
- Suitable Lid cover
- Polishing Multi Strip Holding Vice.
- Cold-finger condensor for IP 227
- Glass cradle for IP 227.



Corossion Accessory Kit available

Country:

Distributor:

