# Flash and Fire Points by Cleveland Open Cup

#### test method

For flash and fire points of all petroleum products, except fuel oils and those having an open cup flash below 79°C (175°F).

### Cleveland open cup flash tester

- · Conforms to ASTM D92 and related specifications
- For flash points above 79°C (175°F)

Determines flash and fire points by the Cleveland Open-Cup method. Consists of test flame applicator, brass test cup, thermometer support, heating plate and electric heater. Applicator is precisely aligned per specifications and pivots for test flame application at specified temperature intervals. Hinged thermometer support raises to facilitate placement and removal of test cup. Adjust flame size using built-in needle valve and comparison bead.

Equipped with a 1000W nickel-chromium heater with stepless variable heat control for accurate repeat setting of temperature rate of rise per specifications.



#### ordering information

description catalog no.

K13902 Electrically Heated CE Model, 115V 60Hz K13992 Electrically Heated CE Model, 220-240V 50/60Hz

accessories

**250-000-11F** ASTM 11F Thermometer Range: 20 to 760°F 250-000-11C ASTM 11C Thermometer Range: -6 to +400°C K14000

Cleveland Open Flash Cup

Precision machined brass with heat resistant handle

## specifications

Conforms to the specifications of:

ASTM D92; D6074, D6158; AASHTO T48; ANS Z- 11.6; IP 36; ISO 2592; DIN 51376; FTM 791-1103, FTM 141-4294

**Electrical Requirements:** 

115V 60Hz

220-240V 50/60Hz, Single Phase, 3.4A

Dimensions wxdxh in.(cm)

10 x 17.5 x 11.5 (25.4 x 44.5 x 29.2)

Net Weight: 20 lbs (9.1kg)

**Included Accessories** 

Brass Test Cup

**Shipping Information** 

Shipping Weight: 30 lbs (13.6 kg) Shipping Dimensions: 2 Cu. ft.

